


GE Oil & Gas – Reuter-Stokes

GE Reuter-Stokes

Jonathan Voss
Dr. Mathieu Boucher

June 13, 2017



Our People - Twinsburg, OH Facility



OSHA VPP Star
GE eCO2 certified
HealthAhead
ISO 9001
10 CFR 50 App B
ASME Sec III
200 EEs
Non-union
163'000 ft²
100 acres

Reuter-Stokes key milestones

1956 – Reuter-Stokes founded

1984 – Reuter-Stokes acquired by GE

2001 – Tensor (Oil & Gas exploration products) acquired by GE

2008 – ^3_2He shortage

2010 – Commercial ^3_2He available

2010 – Reuter-Stokes launches $^{10}_5\text{B}$ lined proportional counter


2011 – Reuter-Stokes investigates $^{10}_5\text{B}$ for safeguards

Reuter Stokes Measurement Solutions

Technologies for Harsh Environments

Nuclear Instrumentation

Neutron Monitoring
In-core power monitoring
Ex-core power monitoring
Gamma Thermometer



Radiation Measurements

Homeland Security
Safeguards
Scattering
Environmental monitoring
 ^3He recycling




Drilling Tools

^3He Detectors
Oil exploration
Geo-steering



Turbine Instruments

Flame detection sensor
- SiC technologies




Reuter Stokes Measurement Solutions

Technologies for Harsh Environment

Nuclear Instrumentation

Neutron Monitoring
In-core power monitoring
Ex-core power monitoring
Gamma Thermometer



Radiation Measurements

Homeland Security
Safeguards
Scattering
Environmental monitoring
 ^3He recycling




Drilling Tools

^3He Detectors
Oil exploration
Geo-steering



Turbine Instruments

Flame detection sensor
- SiC technologies



Radiation Measurement Solutions

60 Years of Experience

- Over 10,000 neutron counter designs manufactured
- ^3He gas-filled detectors – multiple designs
- ^{10}B lined detectors – reactor power monitoring
- BF_3 gas-filled detectors
- High pressure ion chamber
- Scintillator packaging for harsh environments (O&G)
- ^3He recovery and purification

Core Competencies

- Harsh Environment Packaging
- Handling Sensitive Materials
- Radiation Detection
- Precision Fabrication
- Product Testing & Modeling



RS Skills

- Brazing, welding
- Long lived
- Integrated electronics
- Radioactive Sources
- Hygroscopic materials
- NRC compliance / DOT shipping regulations
- neutron
- γ radiation
- Ultraviolet
- Thin metal welding
- Coatings – electroplating and boron
- thin wire and cable processing
- Environmental testing
- Position calibration
- Monte Carlo simulation



Radiation Measurement Solutions

60 Years of Experience

- Over 10,000 neutron counter designs manufactured
- ^3He gas-filled detectors – multiple designs
- ^{10}B lined detectors – reactor power monitoring
- ~~BF_3 gas-filled detectors~~
- High pressure ion chamber
- Scintillator packaging for harsh environments (O&G)
- ^3He recovery and purification




Core Competencies

- Harsh Environment Packaging
- Handling Sensitive Materials
- Radiation Detection
- Precision Fabrication
- Product Testing & Modeling

RS Skills

Failure rate for Portal Detectors




MTBF: 832 years

- Monte Carlo simulation



imagination at work

^3_2He Detector




Detector: Construction


- Metallic outer housing
- Small diameter anode
- Proportional fill gas
- Welded and brazed construction
- Optimized ^3He pressure for desired efficiency
- Optimized size for instrument applications

Operation:

- Well understood principles
- Neutron captured by ^3He nucleus
- Gamma response primarily from interaction with detector wall
- Charge sensitive electronics
- Well defined spectrum



^3_2He Detector




Detector: Construction

- Metallic outer housing
- Small diameter anode
- Proportional fill gas
- Welded and brazed construction
- Optimized ^3He pressure for desired efficiency
- Optimized size for instrument applications

Operation:

- Well understood principles
- Neutron captured by ^3He nucleus
- Gamma response primarily from interaction with detector wall
- Charge sensitive electronics
- Well defined spectrum



GE Measurement & Control – ^3_2He Detectors

Available Sizes

Detector Diameter:

- ↳ from 8mm
- ↳ up to 5 inches (~ 11cm)

Detector sizes (active length):

- ↳ 8 mm → up to 41 inches
- ↳ $\frac{1}{2}$ inch → up to 44 inches
- ↳ 1 inch → up to 150 inches

Aluminum: $\min \varnothing = \frac{1}{2}''$

Fill pressure

0.1 atm \leq **stainless steel** \leq 30 atm

The maximum fill pressure for **aluminum** is 15 bar.

For shipping purposes, all detectors must be certified up to 4 × the fill pressure.

Materials of construction

Stainless steel:

- ↳ Background: $2.9 \times 10^{-4} \frac{\text{cps}}{\text{cm}^2}$

Aluminum:

- ↳ Background: $2.9 \times 10^{-3} \frac{\text{cps}}{\text{cm}^2}$

Detector Type

- ↳ Proportional Counter
- ↳ PSD

Connector

- ↳ Vacuum Couplings
- ↳ straight or right angle connector
- ↳ HN, MHV, SHV, Lead wire

Every detector is custom made



NeuAcq® Overview

Introduction and History
System Overview
8-pack Overview
Networking Standards
Communication Protocols



NeuAcq® Overview

Introduction and History

- 2007 GE's Reuter-Stokes licensed electronics from SNS,
- Based on non-proprietary, existing network standards,
 - ↳ achieve higher data rates,
 - ↳ improved reliability,
 - ↳ significant simplification in the overall system interconnections and cables, and
 - ↳ off the shelf network devices.

System Overview

8-pack Overview

Networking Standards


Communication Protocols

NeuAcq® Overview

Introduction and History

System Overview

- Fully scalable system and electronics
 - ↳ blue components ⇒ off-the-shelf
 - ↳ mitigated obsolescence
- Network Routers
 - ↳ IEEE 1588-2008 PTP,
 - ↳ IEEE 802.3af POE,
 - ↳ No. of routers depends on system size.
- Power ⇒ PoE
- Grandmaster clock
 - ↳ time base for the PTP functionality,
 - ↳ typically uses a GPS clock.
- Timestamper
 - ↳ provided by GE,
 - ↳ external synchronization for ToF,
 - ↳ absolute time measurement,
 - ↳ version of the RSPP



NeuAcq® Overview

Introduction and History System Overview

8-pack Overview

- ^3He PSD,
 - ↳ charge division,
 - ↳ various detector lengths, diameter, and sensitivity.
- Decoupling Capacitors,
 - ↳ HV isolation, while passing pulses
- Preamplifiers,
- ADCs,
 - ↳ very fast and highly accurate,
- ↳ digitizing pulses at 10^7 samples/s.
- Platform Processor,
 - ↳ RSPP \Rightarrow FPGA and ARM processor.
- Waveform Digitization can compensate:
 - ↳ baseline drift,
 - ↳ detection of multiple pulses, and
 - ↳ γ /neutron discrimination.

Networking Standards Communication Protocols

NeuAcq® Overview

Introduction and History

System Overview

8-pack Overview

Networking Standards

- POE,
 - ↳ IEEE-802.3.at-2009 (POE+)
 - ↳ up to 25.5 W for 48-44 VDC range,
 - ↳ no need for extra power cable,
 - ↳ power isolated \Rightarrow minimize ground loop issues
- PTP,
 - ↳ very accurate and stable timestamp over Ethernet,
 - ↳ GPS based grandmaster clock,
 - ↳ absolute time stamping.
- TCP,
 - ↳ all event data,
 - ↳ guaranteed correct data delivery,
 - ↳ tests with 1GB Ethernet
 - $\Rightarrow 6 \times 10^6 \frac{\text{events}}{\text{second}}$,
 - ↳ higher rates could be achieved with faster network.
- UDP.
 - ↳ discover mechanism \Rightarrow which modules are on network,
 - ↳ could be eliminated with static IP address.

Communication Protocols

NeuAcq® Overview

Introduction and History

System Overview

8-pack Overview

Networking Standards

Communication Protocols

Compatible with:


- EPICS
- ADARA



GE's Reuter-Stokes partnership with JCNS-KWS-2

- Development of Multi-MHz detector array, in partnership with JCNS
 - ↳ 18 × 8-pack of three (3) different length,
 - ↳ 2 × GB routers,
 - ↳ Grandmaster clock,
 - ↳ timestamplers, and
 - ↳ mounting frame

<http://www.fz-juelich.de/SharedDocs/Meldungen/JCNS/EN/2015/2015-11-06-SANS-2-Detector.html>




imagination at work

GE Reuter-Stokes
13 / 20
June 13, 2017

GE's Reuter-Stokes partnership with JCNS-KWS-2

– Development of Multi-MHz detector array, in partnership with JCNS

- ↳ 18 × 8-pack of three (3) different length,
- ↳ 2 × GB routers,
- ↳ Grandmaster clock,
- ↳ timestamps, and
- ↳ mounting frame



<http://www.fz-juelich.de/SharedDocs/Meldungen/JCNS/EN/2015/2015-11-06-SANS-2-Detector.html>



imagination at work


GE's Reuter-Stokes partnership with JCNS-KWS-2

System Characteristics

- 0.9m² array of 8mm ^3He detectors
 - ↳ $\eta \sim 85\%$ for 5Å
- dead-time constant of 25ns
- count rate as high as 5MHz with 10% dead-time
- 25× improvement compared to old detector

To improve the read-out characteristics and reduce noise, the electronics is mounted in a closed container at the rear of the detector.

<http://www.fz-juelich.de/SharedDocs/Meldungen/JCNS/EN/2015/2015-11-06-SANS-2-Detector.html>



GE's Reuter-Stokes partnership with JCMS-KWS-2

System



- 0.9m
- ↳ η
- dead
- count
- time
- 25x

To
red
clo

GE Reuter-Stokes measurement solutions

- individual He-3 detectors
- individual B-10 detectors
- individual B-10+ detectors
- detector array, with mounting frame
- NeuAcq® data acquisition electronics

GE's Reuter-Stokes can assist for your measurement solution, from individual detectors to turnkey solutions.




<http://www>
2015-11-06



imagination at work

$^{10}_5\text{B}$ Detector



$$\sigma = 3840\text{ barns @ } 0.025\text{ eV}$$

Detector: Construction similar to ^3He

- Metallic outer housing
- Small diameter anode
- Proportional fill gas
- Welded and brazed construction

$^{10}_5\text{B}$ Lining: organic formulation

- Thin layer on internal shell
- Optimize coating thickness
- Use boron enriched in $^{10}_5\text{B}$ isotope

Operation: new design

- Same physics principles as ^3_2He
- Charge sensitive electronics
- Well defined spectrum
- Validated through use in nuclear reactor instrumentation




$^{10}_5\text{B}$ Detector

$^{10}_5\text{B}$ Coating Optimization

Thickness: optimize neutron counting

Too thin: fewer neutron reactions/higher probability of reaction products escaping into detector volume

Too thick: more neutron reactions/reaction products absorbed in lining





B10Plus+® Neutron Detector

How does it work?

B10Plus+® is :

GERs $^{10}_5\text{B}$ lined proportional counter
incorporates a small amount of ^3_2He = boost neutron sensitivity






B10Plus+ “Superpack”

B-10 lined detector + $\frac{1}{2}$ atm ^3He

4 rows deep

Optimized design using GEANT4

- ↳ Coating thickness
- ↳ Gas fill composition
- ↳ Inter-column shielding



45.4% absolute efficiency (PuBe source, 25 mm collimator) equivalent to
3 atm He-3 8-pack

GERS outlook

Gamma Scintillator Packages for downhole operation

CeBr₃ crystals

NeuAcq® upgrade

- ↳ higher count rate capability
- ↳ more channels per platform processor ⇒ cheaper system

change in GE philosophy



- ↳ possibility to lease larger system from GE's Reuter-Stokes

GE's Reuter-Stokes recognizes the need for better position resolution in scattering detector arrays

Neutron Coincidence Counter

B10Plus+® UNCL

- Four footprints were considered,
- Cavity dimensions configured for PWR fuel rod bundles ($23.4\text{cm} \times 23.4\text{cm}$).
- 85 configurations were simulated using GEANT4, with a ^{252}Cf source at CoG,
- The number of detectors and quantity of ^3_2He gas were minimized to reduce cost and conserve ^3_2He .




(a) existing UNCL design.

(b) based on the UNCL, with up to 2 additional rows of detectors.

(c) and (d) are designed to maximize the number of detectors nearest the cavity.

Variables in GEANT4 simulations

- Detector diameter
 - ↳ 25 mm – 28.6 mm
 - ↳ **25mm**
- Moderator thickness
 - ↳ 9 mm – 14 mm
 - ↳ **14mm**
- Thickness of moderator in front of first row
 - ↳ 6.35 mm – 19.1 mm
 - ↳ **12.7mm**
- Thickness of moderator between rows
 - ↳ 6.35 mm – 12.7 mm
 - ↳ **12.7mm**
- Active length
 - ↳ 14 inch – 24 inch
 - ↳ **18 inch**
- ^3He partial pressure
 - ↳ 0 atm – 1 atm
 - ↳ 0.5 atm (inner row), 0.25 atm (outer row), 4.9 litres total
- ^{10}B coating loading
 - ↳ 0.1 mg/cm² – 0.55 mg/cm²
 - ↳ **0.24 mg/cm²**



UNCL design.

on the UNCL, with additional rows of de-

are designed to maximize number of detectors cavity.



GE imagination at work